

E H M E D I C A L

Product catalogue summary



Keep Moving *Forward*™



Featuring
Motio StepWatch™
and
Motio Functional Report
technology

Upper extremity:

- Essential Shoulder Immobilizer
- Deluxe Shoulder Immobilizer
- Pediatric Shoulder Sling
- Apollo Universal Wrist brace 8" & 10"
- Pediatric Universal Wrist brace
- Wrist Lacer with/out Thumb spica




Lower extremity:

- T Scope Premier Post op Knee brace
- Breg Ankle Wraptor
- Genesis Walkers
- PFS Soft Night Splint
- PFS Straps
- Deluxe- and Square toe post op shoe
- Square Toe Offloading post op shoe



Spinal braces:

- Basic Lumbar support




Functional Braces:

Axiom-D Elite Ligament Knee Brace

The Axiom-D Elite's unique dynamic hinge delivers a progressive counterforce to the lower leg. This reduce tibial translation in patients with ACL & PCL ligament deficiencies or recovering from ligament reconstruction.

This maintains the proper anatomic relationship between the femur and tibia as it goes through a full range of motion.



DUO Knee Brace

The DUO (Dynamic Unloading Osteoarthritis) brace is the only dual-upright dynamic OA brace. This brace creates a load across the joint when the knee is extended and reduces the load as the knee flexes. Thus, when the patient is bearing weight, the brace is actively unloading. DUO reduces force when sitting, making it ideal for patients who need to wear a brace for daily activities.



Z13 Knee Ligament Brace

The lightweight brace ideal for patients of shorter stature. Indicated in knee ligament injuries / deficiencies and or in Meniscal instabilities





Cold Therapy

Polar Care Kodiak

One of the best-selling cold therapy products today, designed for convenience and versatility.

The whisper-quiet pump technology allows patients to rest and recover in comfort. To be used with the revolutionary Intelli-Flo pads that offer temperature specific to each treatment area.



Polar Care Wave

Combines two therapies into one.

The motorized **cold therapy** with active compression in one simple and compact system.

Adjustable **compression therapy**:

- 50 mmHg or
- 25 mmHg average peak compression

Adjustable cold therapy settings:
7°C or 10°C average operating temperature.



Motio StepWatch™ and Motio Functional Report technology

The problem statement:

- There is a need to quantify the activity level and potential benefit of a prosthetic device for lower-limb amputees.
- Traditional tests can be subjective, and ratings are often based on short, once-off assessment.
- It is difficult to determine the optimal componentry for each prosthetic user needs.
- Reimbursement data are created manually.

The Solution:

Giving funders the real-world data they require to support your patients.

Motio StepWatch™ device



Accompanies a prosthetic user for a minimum period of 7 days, gathering data on cadence, distance, speed, and other metrics to allow the clinician to properly classify their patient's activity and K-level.

Motio Functional Report




Allows O&P professionals to perform current standard outcome measures while assessing the patient's ambulation data acquired from their real life and outside the clinical environment.



Accurate Patient Assessment

To centralize all the essential information required for an accurate patient assessment, the Motio System incorporates different outcome measures that can be used for reimbursement purposes & justifying a patient's need for certain prosthetic components.





E H M E D I C A L

CONTACT US:

Elrise Holtzhausen

+27 72 630 9911

elrise@ehmedical.co.za